



LINKMED DICOM
Web-enabled Interface Engine

HL7, XML to DICOM and DICOM to HL7, XML *Conversion*

LINKMED DICOM Features:

- HL7 TCP/IP Receiver
- HL7 TCP/IP Transmitter
- Dynamic Message Mapper
- Web-enabled Interface Engine
- Interface Maintaining Utilities

LINKMED DICOM provides the following interfaces:

HL7

- Receives HL7 ADT, Order, Scheduling messages from HIS/EMR systems
- Converts HL7 messages to DICOM files
- Converts DICOM Files to HL7 messages
- Sends HL7 ORU messages to HIS/EMR systems

XML

- Converts XML to DICOM file
- Converts DICOM file to XML

System Requirements:

Windows® 2000/XP Pro/Vista or
Windows® Server 2000/2003 or above
Pentium 4, AMD 64 2.0 GHZ or above CPU
100 GB Hard Disk Space
1 to 2 GB of RAM or more